

# The perfect contact lens for exceptional vision!



*Vision uncompromised.*

**Duette<sup>®</sup> HD**  
CONTACT LENSES

# Why should your patients settle for an average lens?

- **Almost half** of toric wearers report that their vision is blurry or hazy
- **Over a third** say they have fluctuating vision
- **A third** have experienced glare or halos in low light conditions



With Duette HD lenses, even night driving is crystal clear.

## Are those patients looking for better vision elsewhere?

### Duette patients are happy patients!

#### Uncompromised Optics – Even at Night

- Aspheric GP center delivers exceptional vision, even in low illumination, with greater depth of focus

#### Stable Vision Throughout the Day


- Vision not affected by lens rotation

#### Excellent Ocular Health


- Hyper Dk materials provide excellent oxygen permeability
- UV blocker protects from the sun's harmful rays

#### Comfort by Design

- All-day tear exchange and lens movement make Duette lenses extremely comfortable



#### How do Duette HD prescribers rate the lens?



*More than 140 Duette HD prescribers were asked to rate visual stability and visual acuity on a scale of 1 to 5.*

## A new premium option for your contact lens patients.

You have many premium options for your spectacle patients seeking better vision, but your patients wearing contacts haven't had a superior alternative to commodity products...until now.

Duette HD is the ideal lens for patients who want the best possible vision; vision that average contact lenses do not deliver. In fact, 92% of prescribers surveyed rated the vision of Duette HD superior to the vision of soft torics.

**Athletes and weekend warriors demand a lot from their vision.**



## Patient candidates – your everyday patients.

**Duette HD is a great lens for any patient who wants the very best acuity.**

#### Reach for Duette HD for your astigmatic patient who:

- Has decreased or fluctuating vision
- Is wearing glasses because contacts did not deliver acceptable vision
- Has mild to moderate corneal astigmatism or corneal cyl less than 0.75D
  - Is an emerging presbyope with an add of 1.00D or less
  - Would like better vision at night and in low light conditions

*Vision uncompromised.*

**Duette<sup>®</sup> HD**  
CONTACT LENSES



## For independent eye care professionals only: Exclusive practice benefits

We're committed to partnering with you by bringing you products that are available exclusively through the independent ECP. This gives you an opportunity to differentiate your practice with a premium contact lens that delivers exceptional vision. And since Duette HD contact lenses are not available online or through mass merchandisers, you'll build a profitable annuity for your practice, as patients return year after year.

### Fitting Duette HD is Easy!

Duette HD lenses may be fit empirically without the need for diagnostic sets or fluorescein. And the new alignment fitting approach produces consistent, predictable results and successful outcomes.



**Duette HD is ideal for people  
who need to see their best.**

Duette HD	Parameters
Diameter	14.50mm
Base Curve	7.10 to 8.30 in 0.20 steps
Skirt Curve	8.7 (Flat2), 8.4 (Flat), 8.1 (Medium)
Sphere Power	+10.00 to -15.00 <ul style="list-style-type: none"><li>• +8.00 to -8.00 in 0.25D steps</li><li>• +8.50 to +10.00 in 0.50D steps</li><li>• -8.50 to -15.00 in 0.50D steps</li></ul>
Materials	130 Dk GP center, 84 Dk SiHy skirt
Class II UV blocker	>80% UVA and 95% UVB radiation

### Empirical Ordering

**USA & Canada | Monday-Friday: 9:00AM to 8:00PM (EST)**

#### Technical Consultation

Tel: 877.733.2012, option 2 E-Mail: [consultation@synergieyes.com](mailto:consultation@synergieyes.com)

#### Place Orders & Customer Care

Tel: 877.733.2012, option 1

Fax: 877.329.2012

E-Mail: [orders@synergieyes.com](mailto:orders@synergieyes.com)

Online: [www.synergieyes.com](http://www.synergieyes.com)

**Outside USA & Canada |** To locate an International Distributor visit [www.synergieyes.com/international](http://www.synergieyes.com/international).

To contact Customer Care in USA call +1 760.476.9410, option 1, or to place an order email: [intorders@synergieyes.com](mailto:intorders@synergieyes.com).

An Empirical calculator is also available on [SynergEyes.com](http://SynergEyes.com) in the lens ordering area.