



UVEOSCLERAL IMPLANT V-CLIP



Product's Indications Description

•The **ESNOOPER CLIP** is an Uveoscleral implant specifically designed for the technique Non-perforating **Deep Sclerectomy** in open angle glaucoma.

• With the new implant **ESNOOPER CLIP** we claim to take advantage and promote two differentiated drainage tracts, **TRABECULAR** (INTRASCLERAL and SUPRACILIARY) and **UVEOSCLERAL**, inserting a plate in the suprachoroidal bag and positioning the other over the scleral bed, flowing the aqueous humour through them.

PRODUCT'S SPECIFICATIONS

Comercial Name	ESNOOPER CLIP
Product's code	VCLIP-AJL
Definition	Uveoscleral implant: Esnoper
Material	Hydroxyethyl methacrylate (HEMA)
Sterilization method	Steam
Expiration Date	3 years



Presentation

• The implant **ESNOOPER CLIP** is in sterile single packs, easy to open, immersed in highly purified water, and placed in special forceps so that the implant is visible and has an optimum access for its handling,

Quality Guarantee

• Meets the requirements of **Medical Devices Directive 93/42/CEE** and surpassed to the national legislation in the decree **RD 1591/2009** and the Standards **ISO 9001, ISO 13485** "Medical Devices. Quality Management System".

• **AJL** holds the certifications **ISO 13485, ISO 14001** and **ISO 9001**

PRODUCTION AND DISTRIBUTION: www.ajlsa.com

